Automatic Blood Pressure Monitor Model M1-2

INSTRUCTION BOOK

Thank you for choosing our Automatic Blood Pressure Monitor Gamma™ designed on the basis of the most advanced technologies. The device you are holding now is your reliable assistant. You can measure your blood pressure with this monitor as accurately and precisely as an experienced doctor or a skilled nurse would do it.

We are sure that having appraised worthily the high quality and reliability of this device you will become a regular user of the products of Gamma Trade Mark.

Before starting to use this device please study the instruction book carefully. The instruction book offers all information you need to correctly measure your blood pressure and pulse.

Regular blood pressure control allows for reliable and easy health monitoring and for effective combating possible complications caused by blood pressure abnormalities.

ATTENTION

This blood pressure monitor is designed to carry out self-control over blood pressure but NOT to make self-diagnosis of hypertension/hypotension. Please DO NOT diagnose by yourselves basing on the measurement results obtained with the blood pressure monitor. Please DO NOT execute self-treatment of high/low blood pressure and DO NOT change the methods prescribed without consulting your doctor.

CONTENT

1. Features of Automatic Blood Pressure Monitor Model M1-2

- 1.1. Components of your blood pressure monitor kit
- 1.2. Description of the blood pressure monitor

2. Getting ready for measurement

- 2.1. Inserting the batteries
 - 2.1.1. Using an A/C power adapter
- 2.2. Cuff connection
- 2.3. Setting the time and date

3. Measurement Procedure

- 3.1. Fitting the cuff
- 3.2. Measuring procedure
- 3.3. Reading measurement results
- 3.4. Discontinuing a measurement
- 3.5. Irregular Heartbeat Detector

4. Memory function

5. Error Messages/Troubleshooting

6. Care and maintenance

7. Technical specifications

8. Warrantee

9. Important information on blood pressure and its measurement

- 9.1. What is blood pressure?
- 9.2. Which values are normal?
- 9.3. What can be done to lower/maintain the level of blood pressure?


1. Features of Automatic Blood Pressure Monitor Model M1-2

1.1. Components of your blood pressure monitor kit

The blood pressure monitor kit includes Automatic Blood Pressure Monitor Model M1-2, the cuff, 4 batteries (AA), Instruction Book, Warranty Card, kit box.

1.2. Description of the blood pressure monitor

Measuring unit Model M1-2


2. Getting ready for measurement

2.1. Inserting the batteries

To insert batteries

1. Remove the cover as illustrated.

2. Insert the batteries (4 AA, 1.5 V), following the indicated polarity. (+ -)

Attention

- If the battery warning  appears in the display, the batteries are almost flat and must be replaced by new ones
- After the battery warning  appears, the device is blocked until the batteries have been replaced.
- If the blood pressure monitor is not used for long periods, remove the batteries from the device.
- Please use "AA" Long-Life or Alkaline 1.5V Batteries.

Do not use rechargeable batteries.

- Functional check: Press M button down to test all the display elements. When functioning correctly all segments must appear

2.1.1. Using an A/C power adaptor


It is possible to operate Automatic Blood Pressure Monitor Model M1-2 with an AC/DC adaptor.

1. Push the plug into the socket at the bottom part of the device. No power is taken from the batteries while the AC/DC adaptor is connected to the device.

2. Plug the AC adaptor into a 220/110 V power socket. The de-

vice is ready for measurement.

3. Test that power is available by pressing the START button.


Note:

- Please use the AC adaptor of trade mark Gamma. Any other adaptor being used, the service centre is not considered responsible for the warranty servicing/repair of the monitor.
- No power is taken from the batteries while the AC/DC adaptor is connected to the device.
- Please consult if you have questions relating to the AC/DC adaptor.

2.2. Cuff connection

Insert the cuff tube into the opening provided on the side of the device as shown in the picture.


2.3. Setting the time and date

This blood pressure monitor incorporates the function of time and date display. Time and date are recorded along with pressure and pulse.

After new batteries have been inserted, you come into direct time (Year) setting from: 2005-01-01 00:00 O'clock. For this, please proceed as follows:

1. When new batteries have been inserted, the display comes into Year setting, during which the 4 characters blink.

2. Enter the correct year by pressing the MEMORY button. Each click will change the clock 1 year.

3. Press the TIME button. The display now switches to the current month.

4. Enter the current month by pressing the MEMORY button.

5. Follow the same steps to set the current day.

6. Press the TIME button again. The display now switches to the current time, in the format Hour-Minutes.

7. Enter the corresponding time by pressing the MEMORY button.

8. After all settings have been made, click the TIME button once. The date is briefly displayed and then the time. The input is now confirmed and the clock begins to run.

Note

- If you want to change the date and time, press and hold the time button down for approx. 3 seconds till Year flashes, now you can enter the new values as described above.
- 1 x press on the MEMORY button and the TIME button advances one operation (e.g. switch to the next setting mode or change the figure into +1)
- Holding the button down speeds up the procedure.


3. Measurement Procedure

ATTENTION

- Find time to relax by sitting in a quiet atmosphere for some time before measurement.
- Efforts by the patient to support the arm can increase the blood pressure. Make sure you are in a comfortable, relaxed position and do not activate any muscles in the arm during measurement.
- Always measure on the same arm (normally left).
- Remove any garment that fits closely to your upper arm. Do not roll the sleeve since it can squeeze your hand and this can lead to false results.
- Use only clinically approved original cuff.
- If you want to follow the results of your blood pressure measurements, always perform measurements at the same time of day, since blood pressure changes during the course of day.
- Measurements should be done after a 5 minute rest to ensure accuracy.

3.1. Fitting the cuff

1. Pass the end of the cuff through the flat metal ring so that a loop is formed, and the velcro faces outwards.

2. Place the cuff over the left upper arm

so that tube is closer to your lower arm.

3. Lay the cuff on the arm so that the lower edge of the cuff lies approximately 2 to 3 cm above the elbow.

4. Tighten the cuff by pulling the end and close the cuff by affixing the velcro. There should be little free space between the arm and the cuff. Cuff that does not fit properly results in false measurement values. Clothing must not restrict the arm. Any piece of clothing which does must be removed.

5. Lay your arm on a table so the cuff is at the same height as your heart. If the upper arm artery lies considerably lower or higher than the heart, a false higher or lower blood pressure will be measured! A variation of 15 cm between cuff and heart level can result in a reading error of + or - 10mm Hg.

Make sure the tube is not kinked.







6. Remain seated quietly for two minutes before you begin the measurement.

3.2. Measuring procedure

To carry out measurement

1. Press the START button. The built-in microprocessor begins to inflate the cuff. In the display, the increasing cuff pressure is continually displayed.

2. After automatically reaching an in-

dividual pressure, the pump stops and the pressure slowly falls. The cuff pressure is displayed during the measurement. When the device has detected your pulse, the heart symbol in the display begins to blink and a beep tone is audible for every pulse beat.

3.3. Reading measurement results

When the measurement has been concluded, a long beep tone sounds. The measured systolic and diastolic blood pressure values, as well as the pulse are now displayed.



The appearance of this symbol ♥ signifies that an irregular heartbeat was detected. This indicator is only a caution. It is important that you be relaxed, remain still and do not talk during measurements.

Note: We recommend contacting your physician if you see this indicator frequently.

To prolong the batteries' life the device switches off automatically if no button is pressed for 1 minute. Otherwise you can switch it off by pressing the START button.

3.4. Discontinuing a measurement

If it is necessary to interrupt a blood pressure measurement for any reason (e.g the patient feels unwell), the START button can be pressed at any time. The device then immediately lowers the cuff pressure automatically.

3.5. Irregular Heartbeat Detector

This function indicates allorhythmic heartbeating. If the symbol of IHD ♥ appears on the display that means that certain abnormality in heart beat frequency was detected during the measurement. In this case, the result may deviate from your normal blood pressure – repeat the measurement. In most cases, this is no cause for concern. However, if the symbol ♥ appears on a regular basis (e.g. several times a week with measurements taken daily) we advise you to tell your doctor. Please show your doctor the following explanation:

Information for the doctor on frequent appearance of the Irregular Heartbeat Detector

This instrument is an oscillometric blood pressure monitor that also analyses pulse frequency during measurement. The instrument is clinically tested.


The symbol ♥ is displayed after the measurement, if pulse irregularities occur during measurement. If the symbol ♥ appears more frequently (e.g. several times per week on measurements performed daily) we recommend the patient to seek medical advice.

The instrument does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

4. Memory function

Automatic Blood Pressure Monitor Gamma™ Model M1-2 automatically stores 99 measured results together with the date and the time when they are obtained.

To see the results of the last measurement press and hold the MEMORY button for at least 3 seconds. Each time you press the memory button an earlier measurement will


be displayed along with time and date (with symbols M2-M99).

If an irregular heartbeat interval was detected in any measurements stored in memory, the icon ♥ will be displayed.

In order to delete all stored readings, depress the MEMORY button for at least 7 seconds, the display will show the symbol CL and 3 short beep sounds will be heard to indicate deletion of stored readings.

5. Error Messages/Troubleshooting

If an error occurs during a measurement, the measurement is discontinued and a corresponding error code is displayed (picture: error no. 2)

Error No.	Possible cause(s) / Solutions
ERR 1	The systolic pressure was determined, but then the pressure in the cuff went down to less than 20 mmHg. Pulse could not be detected. The tube may have loosened after the systolic pressure has been determined.
ERR 2	Unnatural pressure impulses prevent from accurate measuring the blood pressure. Reason: The arm was moved during the measurement. Repeat measurement following the above rules.
ERR 3	If inflation of the cuff takes too long, the cuff is not correctly seated or the hose connection's not tight.

ERR 5	The difference between systolic and diastolic is excessive. Measure again carefully following proper procedures. Consult your doctor if you still get unusual results of measurements.
HI	The pressure in the cuff is too high (over 300 mmHg) OR the pulse is too high (over 200 beats per minute). Relax for 5 minutes and repeat the measurement.*
LO	The pulse is too low (less than 40 beats per minute). Repeat the measurement.*

*Please consult your doctor, if this or any other problem occurs repeatedly

Other possible errors and their solutions

If problems occur when using the device, the following points should be checked:

Malfunction	Remedy
The display remains blank when the device is switched on although the batteries are in place.	1. Check battery installation/ polarity. 2. Remove the batteries and if the display is unusual, then exchange them for new ones.
The pressure does not rise although the pump is running.	Check the connection of the cuff tube and connect properly.

<p>The device frequently fails to measure, or the values measured are too low or high.</p>	<ol style="list-style-type: none"> 1. Fit the cuff correctly on the arm. 2. Before starting measurement make sure that clothing is not exerting pressure on the arm. Take articles of clothing off if necessary. <p>Measure blood pressure again in complete peace and quiet.</p>
<p>Every measurement results in different values, although the device functions normally and normal values are displayed.</p>	<ol style="list-style-type: none"> 1. Please re-read the points listed under "Measuring Procedure" 2. Blood pressure changes constantly. The observed readings may accurately reflect your pressure-minor differences in the readings are not abnormal.
<p>Blood pressure values differ from those measured at the doctor's</p>	<p>Record daily values and consult your doctor.</p> <p>Pressure readings in your doctor's office may be higher due to anxiety.</p>
<p>After the instrument has inflated the cuff the pressure falls very slowly, or not at all.</p>	<p>Check cuff connections.</p>

Note




Blood pressure is subject to fluctuations even in healthy people. Please remember that **comparable blood pressure measurements always require the same time and the same conditions! These are normally quiet conditions.** If you follow the correct procedures described above and still get the fluctuations of blood pressure of more than 15 mmHg and/or you repeatedly hear irregular pulse tones please consult your doctor.

ATTENTION

Should any technical malfunctions arise in the blood pressure monitor, please contact the point of sale where you bought the device or Service Centre Gamma-Medica. **Never attempt to repair the instrument yourself!**

Any unauthorized opening of the instrument invalidates all warranty claims.

6. Care and maintenance

<p>Do not expose the device to either extreme temperatures, humidity, dust or direct sunlight.</p>	
<p>Handle the cuff carefully and avoid all types of stress through twisting or buckling in order not to damage the sensitive air-tight bubble.</p>	
<p>Clean the device with a soft, dry cloth. Do not use gas, thinners or similar solvents. Spots on the cuff can be removed carefully with a damp cloth and soapsuds. Do not submerge the cuff in water!</p>	
<p>Handle the tube carefully. Keep the hose away from sharp edges.</p>	
<p>Do not drop the monitor or treat it roughly in any way.</p>	
<p>Never open the monitor! Otherwise the manufacturer's calibration will be spoiled.</p>	

Note:

According to international standards, your monitor should be checked for calibration every 2 years.

7. Technical specifications

Weight	483 g (with batteries)
Size	124(W) x 205(L) x 81(H) mm
Storage temperature	-5 to +50°C (23°F - 122°F)
Humidity	15 to 85% relative humidity maximum
Operation temperature	10 to 40°C (50°F - 104°F)
Display	LCD-Display (Liquid Crystal Display)
Measuring method	oscillometric
Pressure sensor	capacitive
Measuring range:	
- SYS/DIA pressure	30 to 280 mmHg
- Pulse	40 to 200 per minute
Cuff pressure display range	0–299 mmHg.
Memory	Automatically stores the last 99 measurements.
Measuring resolution	1 mmHg
Accuracy	
- pressure	within ± 3 mmHg
- pulse	± 5 % of the reading
Power source	4 x AA batteries 1,5V
Accessories	Automatic Blood Pressure Monitor Model M1-2, cuff, 4 AA batteries, 1.5 V, Instruction Book, Warranty Card, packing.

*technical alterations are possible!

8. Warranty

Your Automatic Blood Pressure Monitor Model M1-2 is warranted for 3 years from date of purchase. This warranty includes the monitor and the cuff.

The warranty does not apply to damage caused by improper handling, accidents, not following the operating instructions or self-maintained alterations made to the device.


The warranty is only valid upon presentation of the warranty card which was correctly filled in and sealed.


9. Important information on blood pressure and its measurement

9.1. What is blood pressure?

Blood pressure is the pressure that your blood exerts on the vascular walls. Blood pressure is necessary to provide for constant blood flow inside the body. Thanks to it the cells get oxygen that provides for their normal functioning. The heart performs the function of a “pump”, sending blood to the blood vessels. Each heart beat creates a certain level of the blood pressure.

There are 2 kinds of blood pressure: a systolic (upper) one, which corresponds to the heartbeat pushing blood into the arteries; and a diastolic (lower) one, which means the blood pressure between two heartbeats.

Blood pressure is subject to fluctuations during the day even in healthy people. The fluctuations are influenced by a number of factors – time of day, person’s condition, physical or mental activity, environment, etc.

Day-Night Fluctuations of the Blood Pressure

An increase of blood pressure increases the burden onto the heart, affects blood vessels making their walls thick and less elastic. One of the features of the hypertension is its ability to remain unnoticed for the patient at its early stages. That’s why the self-control of the blood pressure is so important. With the illness progressing, headaches and regular dizziness appear, the sight declines, the functioning of vitalses (encephalon, heart, kidneys, blood vessels) breaks down. Without special treatment the complications of hypertension might be kidney damages, breast-pang, paralytic stroke, aphasia, dementia, heart attack and stroke.

9.2. Which values are normal?

The world standard as for the norms of the blood pressure is the Classification* of the World Health Organization (WHO):

Blood Pressure Category	Systolic Blood Pressure (mmHg)	Systolic Blood Pressure (mmHg)	Measures
Hypotension	< 100	< 60	Doctor’s control

Optimum pressure	100 ... 120	60 ... 80	Self-control
Normal pressure	120 ... 130	80 ... 85	Self-control
High normal pressure	130 ... 140	85 ... 90	Self-control, mode of life correction
Mild hypertension	140 ... 160	90 ... 100	Consultation with doctor
Moderate hypertension	160 ... 180	100 ... 110	Consultation with doctor
Severe hypertension	> 180	> 110	Consultation with doctor

* printed with curtailments

- The diagnosis of hypertension requires from the patient to combine medical treatment prescribed by the doctor and mode of life correction.
- People with normal pressure and high normal pressure are recommended to carry out self-control of their tension in order to timely take measures to decrease the blood pressure level down to the optimal one without using any medications.
- For people more than 50 years old high level of systolic blood pressure (higher than 140 mmHg) is more crucial than diastolic pressure.
- Even with blood pressure being normal, people run the bigger risk of hypertension development with advancing age.

ATTENTION

If you have normal results of blood pressure measured under calm conditions but your results are excessively high when measured under the conditions of physical or mental exhaustion, this might be a sign of so called brittle (that is unstable) hypertension.

If you suspect that, please consult your doctor.

When measured correctly, if diastolic blood pressure is more than 120 mmHg, it is necessary to call the doctor immediately.

9.3. What can be done to lower/maintain the level of blood pressure?

Recommendation		How much the systolic blood pressure can be lowered
<p>- Try to maintain the weight which is considered standard for your age. Decrease overweight if any!</p>		<p>5 – 20 mmHg with every 10 kilos of weight reduction</p>
<p>- Avoid salt overuse (it's recommended to use no more than 6 g per day)!</p>		<p>2 – 8 mmHg</p>
<p>- Avoid eating fat food! Introduce to your diet some food that is rich in dietary fiber (fruits, vegetables)</p>		<p>8 – 14 mmHg</p>
<p>- Refuse from smoking outright! - Use only limited amount of alcohol (beer – 350 g, or wine – 150 g, or vodka – 50 g)! - Restrain the usage of caffeine (coffee)!</p>		<p>2 – 4 mmHg</p>
<p>- After medical examination please regularly go in for sports, with the following restrictions:</p> <ul style="list-style-type: none"> • go in for endurance exercises not strengthening exercises. • do not train till complete exhaustion. • if you have any chronic diseases and/or you are older than 40 please consult your doctor before you start going in for sports. <p>- Spend at least half an hour out in the open air every day.</p>		<p>4 – 9 mmHg</p>

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