## **Semi-automatic Blood Pressure Monitor Model M1-S-1**

### **INSTRUCTION BOOK**

Thank you for choosing our Semi-automatic Blood Pressure Monitor Gamma™ designed on the basis of the most advanced technologies. The device you are holding now is your reliable assistant. You can measure your blood pressure with this monitor as accurately and precisely as an experienced doctor or a skilled nurse would do it.

We are sure that having appraised worthily the high quality and reliability of this device you will become a regular user of the products of Gamma Trade Mark.

Before starting to use this device please study the instruction book carefully. The instruction book offers all information you need to correctly measure your blood pressure and pulse.

Regular blood pressure control allows for reliable and easy health monitoring and for effective combating possible complications caused by blood pressure abnormalities.

### **ATTENTION**

This blood pressure monitor is designed to carry out self-control over blood pressure but NOT to make self-diagnosis of hypertension/hypotension. Please DO NOT diagnose by yourselves basing on the measurement results obtained with the blood pressure monitor. Please DO NOT execute self-treatment of high/low blood pressure and DO NOT change the methods prescribed without consulting your doctor.

## **CONTENT**

### **1. Features of Semi-automatic Blood Pressure Monitor Model M1-S-1**

- 1.1. Components of your blood pressure monitor kit
- 1.2. Description of the blood pressure monitor

### **2. Getting ready for measurement**

- 2.1. Inserting the batteries
- 2.2. Cuff connection

### **3. Measurement Procedure**

- 3.1. Fitting the cuff
- 3.2. Measuring procedure
- 3.3. Reading measurement result
- 3.4. Inflating further
- 3.5. Discontinuing a measurement

### **4. Memory function**

### **5. Error Messages/Troubleshooting**

### **6. Care and maintenance**

### **7. Technical specifications**

### **8. Warrantee**

### **9. Important information on blood pressure and its measurement**

- 9.1. What is blood pressure?
- 9.2. Which values are normal?
- 9.3. What can be done to lower/maintain the level of blood pressure?


# 1. Features of Semi-automatic Blood Pressure Monitor Model M1-S-1

## 1.1. Components of your blood pressure monitor kit

The blood pressure monitor kit includes Semi-automatic Blood Pressure Monitor Model M1-S-1, cuff with hose and pumping bulb, 4 batteries (AA), Instruction Book, Warranty Card, kit box.

## 1.2. Description of the blood pressure monitor

Measuring unit Model M1-S-1


## 2. Getting ready for measurement

### 2.1. Inserting the batteries

To insert batteries


1. Remove the cover as illustrated.

2. Insert the batteries (4 AA, 1.5 V), following the indicated polarity. (+ -)


### ATTENTION

- After low battery indicator appears, the device won't function until the batteries have been replaced.
- If the blood pressure monitor is not used for long periods, remove the batteries from the device.
- Please use "AA" Long-Life or Alkaline 1.5V Batteries. Do not use rechargeable batteries.
- Functional check: Hold the ON/OFF button down to test all the display elements. When functioning correctly many icons will appear.

## 2.2. Cuff connection

Insert the cuff tube into the opening provided on the side of the device as shown in the picture.


## 3. Measurement Procedure

### ATTENTION


- Find time to relax by sitting in a quiet atmosphere for some time before measurement.
- Efforts by the patient to support the arm can increase the blood pressure. Make sure you are in a comfortable, relaxed position and do not activate any muscles in the arm during measurement.
- Always measure on the same arm (normally left).
- Remove any garment that fits closely to your upper arm. Do not roll the sleeve since it can squeeze your hand and this can lead to false results.
- Use only clinically approved original cuff
- If you want to follow the results of your blood pressure measurements, always perform measurements at the same time of day, since blood pressure changes during the course of the day.
- Measurements should be done after a 5 minute rest to ensure accuracy.

### 3.1. Fitting the cuff


1. Pass the end of the cuff through the flat metal ring so that a loop is formed.


2. Place the cuff over the left upper arm so that tube is closer to your lower arm.


3. Lay the cuff on the arm so that the lower edge of the cuff lies approximately 2 to 3 cm above the elbow.

4. Tighten the cuff by pulling the end and close the cuff by affixing the velcro. There should be little free space between the arm and the cuff. Cuff that does not fit properly results in false measurement values. Clothing must not restrict the arm. Any piece of clothing which does must be removed.


5. Lay your arm on a table so the cuff is at the same height as your heart. If the upper arm artery lies considerably lower or higher than the heart, a false higher or lower blood pressure will be measured! A variation of 15 cm between cuff and heart level can result in a reading error of + or - 10mm Hg.

Make sure the tube is not kinked.


6. Remain seated quietly for two minutes before you begin the measurement.


### 3.2. Measuring procedure

To carry out measurement



1. Press ON/OFF button, device is switched on, where you will see many segments display for about 2 secs, then there flashes a downward triangle, which means device is automatically lowering pressure left in the cuff, you can hold down quick release valve to speed up the process.


2. Once cuff pressure is completely released, you will hear 3 short beeps, with a steady "0" shows on


3. Grasp the pumping bulb with your free hand and inflate the cuff, during which an upward triangle keeps flashing.

4. Stop pumping when the upward triangle terminates flashing, device then starts to deflate.


5. Lay down the pumping ball and remain in a quiet and relaxed sitting position. The measurement now proceeds on its own and is taken automatically.

If no measurement takes place and all that appears in the display is a flashing, upwards pointing arrow ▲, pump the ball once again and inflate to a higher pressure (see clause Inflating further given below). When the device detects a pulse, a heart symbol in the display starts to flash and a beep is heard for every heart beat.

### 3.3. Reading measurement results


When the measurement has been concluded, a long beep tone sounds. The measured systolic and diastolic blood pressure values, as well as the pulse are now displayed, and there flashes a downward triangle indicating people to release the rest of air pressure in the cuff.

The cuff pressure can now be quickly released by pressing the quick release valve at the pumping bulb.


### Note

To prolong the batteries' life the device switches off automatically if no button is pressed for 3 minutes. Otherwise you can switch it off by pressing the ON/OFF button.

### 3.4. Inflating further

If the cuff is not inflated sufficiently, the measurement is stopped after a few seconds and a flashing arrow pointing upwards appears in the display. It is then necessary to inflate at least 20 mmHg higher than the previous value till cuff pressure displays again on the LCD. The instruction to inflate further can


appear several times if the pressure is still not sufficient.


**ATTENTION:**

If the cuff is inflated too far (over 300 mmHg), **H<sup>i</sup>** appears in the display and a warning signal is sounded. The cuff pressure must then be immediately released using the quick-release valve at the pumping bulb!


**3.5. Discontinuing a measurement**

If it is necessary to interrupt a blood pressure measurement for any reason (e.g. the patient feels unwell), the quick-release valve at the pumping bulb can be pressed at any time. The device then immediately lowers the cuff pressure automatically.


**4. Memory function**

Semi-automatic Blood Pressure Monitor Gamma™ Model M1-S-1 automatically stores the result measured.


To see the result of the last measurement, press the MEMORY button

## 5. Error Messages/Troubleshooting

If an error occurs during a measurement, the measurement is discontinued and a corresponding error code is displayed (picture: error no. 2)

Error No.	Possible cause(s) / Solutions
ERR 1	The systolic pressure was determined, but then the pressure in the cuff went down to less than 20 mmHg (diastolic blood pressure can't be measured). The tube may have loosened after the systolic pressure has been determined. Further possible reason: Pulse could not be detected.
ERR 2	Unnatural pressure impulses influence the measurement result. Reason: the arm was moved during the Measurement (Artefact)
ERR 5	The measured reading indicated an unacceptable difference between systolic and diastolic pressures. Take another reading following directions carefully. Contact your doctor if you continue to get unusual readings.

### Other possible errors and their solutions

If problems occur when using the device, the following points should be checked:

Malfunction	Remedy
The display remains blank when the device is switched on although the batteries are in place.	<ol style="list-style-type: none"> <li>1. Check battery installation/ polarity.</li> <li>2. Remove the batteries and if the display is unusual, then exchange them for new ones.</li> </ol>
The pressure does not rise although you are inflating the bulb.	Check the connection of the cuff tube and connect properly.
The device frequently fails to measure, or the values measured are too low or high.	<ol style="list-style-type: none"> <li>1. Fit the cuff correctly on the arm.</li> <li>2. Before starting measurement make sure that clothing is not exerting pressure on the arm. Take articles of clothing off if necessary.</li> </ol> <p>Measure blood pressure again in complete peace and quiet.</p>
Every measurement results in different values, although the device functions normally and normal values are displayed.	Please re-read the clause Measuring procedure above.
Blood pressure values differ from those measured at the doctor's.	Record daily values and consult your doctor.

### Note

Blood pressure is subject to fluctuations even in healthy people. Please remember that **comparable blood pressure measurements always require the same time and the same conditions!** These are normally quiet conditions. If you follow the



correct procedures described above and still get the fluctuations of blood pressure of more than 15 mmHg and/or you repeatedly hear irregular pulse tones please consult your doctor.


### ATTENTION

Should any technical malfunctions arise in the blood pressure monitor, please contact the point of sale where you bought the device or Service Centre Gamma-Medica. **Never attempt to repair the instrument yourself!**

Any unauthorized opening of the instrument invalidates all warranty claims.

### 6. Care and maintenance

<p>Do not expose the device to either extreme temperatures, humidity, dust or direct sunlight.</p>	
<p>Handle the cuff carefully and avoid all types of stress through twisting or buckling in order not to damage the sensitive air-tight bubble.</p>	
<p>Clean the device with a soft, dry cloth. Do not use gas, thinners or similar solvents. Spots on the cuff can be removed carefully with a damp cloth and soapsuds. Do not submerge the cuff in water!</p>	
<p>Handle the tube carefully. Keep the hose away from sharp edges.</p>	

Do not drop the monitor or treat it roughly in any way.	
Never open the monitor! Otherwise the manufacturer's calibration will be spoiled.	

**Note:**

According to international standards, your monitor should be checked for calibration every 2 years.

**7. Technical specifications**

Weight	256 g (with batteries)
Size	106(W) x 110(L) x 65(H) mm
Storage temperature	-5 to +50°C (23°F - 122°F)
Humidity	15 to 85% relative humidity maximum
Operation temperature	10 to 40°C (50°F - 104°F)
Display	LCD-Display (Liquid Crystal Display)
Measuring method	oscillometric
Pressure sensor	capacitive
Measuring range:	
- SYS/DIA pressure - Pulse	30 to 280 mmHg 40 to 200 per minute
Cuff pressure display range	0-299 mmHg.
Memory	Automatically stores the result of last measurement.

Measuring resolution	1 mmHg
Accuracy	
- pressure - pulse	within $\pm 3$ mmHg $\pm 5$ % of the reading
Power source	4 x AA batteries 1,5V
Accessories	Semi-automatic Blood Pressure Monitor Model M1-S-1, the cuff with hose and pumping bulb, 4 batteries (AA), Instruction Book, Warranty Card, kit box.

\*technical alterations are possible!

## 8. Warranty

Your Semi-automatic Blood Pressure Monitor Model M1-S-1 is warranted for 3 years from date of purchase. This warranty includes the monitor and the cuff.

The warranty does not apply to damage caused by improper handling, accidents, not following the operating instructions or self-maintained alterations made to the device.


The warrantee is only valid upon presentation of the warrantee card which was correctly filled in and sealed.

## 9. Important information on blood pressure and its measurement

### 9.1. What is blood pressure?


Blood pressure is the pressure that your blood exerts on the vascular walls. Blood pressure is necessary to provide for constant blood flow inside the body. Thanks to it the cells get oxygen that provides for their normal functioning. The heart performs the function of a "pump", sending blood to the blood vessels. Each heart beat creates a certain level of the blood pressure.

There are 2 kinds of blood pressure: a systolic (upper) one, which corresponds to the heartbeat pushing blood into the arteries; and a diastolic (lower) one, which means the blood pressure between two heartbeats.

Blood pressure is subject to fluctuations during the day even in healthy people. The fluctuations are influenced by a number of factors – time of day, person’s condition, physical or mental activity, environment, etc.

### Day-Night Fluctuations of the Blood Pressure

An increase of blood pressure increases the burden onto the heart, affects blood vessels making their walls thick and less elastic. One of the features of the hypertension is its ability to remain unnoticed for the patient at its early stages. That's why the self-control of the blood pressure is so important. With the illness progressing, headaches and regular dizziness appear, the sight declines, the functioning of vitals (encephalon, heart, kidneys, blood vessels) breaks down. Without special treatment the complications of hypertension might be kidney damages, breast-pang, paralytic stroke, aphasia, dementia, heart attack and stroke.

### 9.2. Which values are normal?

The world standard as for the norms of the blood pressure is the Classification\* of the World Health Organization (WHO):

<b>Blood Pressure Category</b>	<b>Systolic Blood Pressure (mmHg)</b>	<b>Systolic Blood Pressure (mmHg)</b>	<b>Measures</b>
Hypotension	< 100	< 60	Doctor's control
Optimum pressure	100 ... 120	60 ... 80	Self-control
Normal pressure	120 ... 130	80 ... 85	Self-control
High normal pressure	130 ... 140	85 ... 90	Self-control, mode of life correction
Mild hypertension	140 ... 160	90 ... 100	Consultation with doctor
Moderate hypertension	160 ... 180	100 ... 110	Consultation with doctor
Severe hypertension	> 180	> 110	Consultation with doctor

\* printed with curtailments

- The diagnosis of hypertension requires from the patient to combine medical treatment prescribed by the doctor and mode of life correction.
- People with normal pressure and high normal pressure are recommended to carry out self-control of their tension in order to timely take measures to decrease the blood pressure level down to the optimal one without using any medications.
- For people more than 50 years old high level of systolic blood pressure (higher than 140 mmHg) is more crucial than diastolic pressure.
- Even with blood pressure being normal, people run the bigger risk of hypertension development with advancing age.

## ATTENTION

If you have normal results of blood pressure measured under calm conditions but your results are excessively high when measured under the conditions of physical or mental exhaustion, this might be a sign of so called brittle (that is unstable) hypertension. If you suspect that, please consult your doctor.

When measured correctly, if diastolic blood pressure is more than 120 mmHg, it is necessary to call the doctor immediately.

### 9.3. What can be done to lower/maintain the level of blood pressure?

Recommendation		How much the systolic blood pressure can be lowered
- Try to maintain the weight which is considered standard for your age. Decrease overweight if any!		5 – 20 mmHg with every 10 kilos of weight reduction
- Avoid salt overuse (it's recommended to use no more than 6 g per day)!		2 – 8 mmHg
- Avoid eating fat food! Introduce to your diet some food that is rich in dietary fiber (fruits, vegetables)		8 – 14 mmHg
- Refuse from smoking outright! - Use only limited amount of alcohol (bier – 350 g, or vine – 150 g, or vodka – 50 g)! - Restrain the usage of caffeine (coffee)!		2 – 4 mmHg

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